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### Instruction Manual Nursing Bed DALI II

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
#### Please notice:

In the text of these instructions square brackets and bold figures = [00] refer to the parts of the nursing bed.

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- 1 Trapeze
- 2 Trapeze pole
- 3 End board, foot section
- 4 Release levers for side rails (4)
- 5 Lower leg rest
- 6 Thigh rest
- 7 Side rail beams (4)
- 8 Hand switch
- 9 Backrest
- 10 End board, head section
- 11 Castors with wheel locks (4)
- 12 Reclining surface frame
- 13 Motor drives for backrest and thigh rest
- 14 Control unit
- 15 Mattress brackets (4)
- 16 Trapeze pole sockets (2)
- **17** Guide bars (4)

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#### **PREFACE**

Dear customer,

Thank you for the confidence you have shown in BURMEIER.

The nursing bed DALI II is factory-checked for electrical safety and function and has left our factory in faultless condition.

Please take time to read this manual carefully.

This applies especially to the assembly phase but also to daily use.

It is also a practical reference guide. Keep the instructions in a handy place.

These instructions provide information for the operator and user regarding the convenient handling and safe operation.

We wish you and your users much success in care-giving and we are convinced that our product will make a positive contribution.

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### **GENERAL NOTES**

These instructions may describe functions and accessories which are not included with your model.

#### Overview DALI II Models:

Model Number: 15.0604.11 Wooden Reclining surface, 4 parts Iron side rails, 2 bars

Model Number: 51.0604.33 Wooden Reclining surface, 4 parts Wooden side rails, 2 bars

Model Number: 15.0604.11 Iron reclining surface, 4 parts Iron side rails, 2 bars

Model Number: 51.0604.33 Iron reclining surface, 4 parts Wooden side rails, 2 bars

### Before putting the nursing bed into operation for the first time:

- Read these instructions carefully in order to avoid damage or incorrect operation.
- Clean and disinfect the nursing bed before first use.

According to the "Medizinprodukte-Betreiberverordnung (MPBetreibV) § 2, before using the nursing bed, the user shall ensure that it is in proper working order and free of defects, and be aware of the instruction manual. This applies also to accessories.

The nursing bed DALI II complies with all requirements of the guidelines 93/42/EWG for medical products. It is classified as a dynamic medical product of the class I in accordance with "Medizinproduktegesetz (German abbreviation: MPG § 13, Medical Products Act)."

The nursing bed DALI II is checked by an independent examining institute. As with all technical, electrical devices improper handling can lead to damage and/or injury.

Observe your obligations as operator in accordance with medical devices "Medizinprodukte-Betreiberverordnung (German abbreviation: MPBetreibV, Operators of Medical Products Ordinance)" in order to ensure a dependable and safe operation of this medical device without endangering patients, users and third parties.

This manual contains safety guidelines, which must be observed. All persons who work with the nursing bed DALI II must familiarize themselves with these instructions and follow the safety guidelines.

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#### **DEFINITIONS OF INVOLVED PERSONS**

In this manual the following persons are involved:

#### Operator

Operators (e.g.: Medical Supply Houses, Dealers and representatives, health insurance companies) are natural persons or legal entities, who utilize the nursing bed DALI II or authorize its use. The operator is responsible for instructing the user in the proper use of the unit.

#### User

Users are persons who as a result of training and experience are authorized to assemble. perform maintenance and/or operate the nursing bed DALI II. Furthermore, the user can recognize and avoid possible dangers, and judge the clinical condition of the patient.

In this manual a patient is an infirm or handicapped person who reclines on the nursing bed.

#### Specialist staff

Specialist staff are employees of the operator who are authorized as a result of their education and training to deliver, assemble, dismantle and transport the nursing bed. Furthermore, they are familiar with the regulations for cleaning and disinfecting the unit.

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#### SAFETY INSTRUCTIONS

The nursing bed DALI II represents the latest technological development at the time of delivery and is checked by an independent examining institute.

■ Do not use the nursing bed DALI II if it is damaged.

#### MEANING OF THE SAFETY SYMBOLS

In these instructions the following safety symbols are used:

#### Warning about injuries to persons

Dangerous voltage. Life Threatening.

General danger. Injury or life hazard.

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#### Warning about property damages

Possibility of damage to drive, material or environment.

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#### Other symbols

Useful tip. For easier operation or better understanding of the unit.

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#### **S**AFETY INSTRUCTIONS FOR OPERATORS

- Remember your responsibilities in accordance with the MPBetreibV in order to ensure the safe operation of this unit without endangering patients, users and third parties.
- Before first use, instruct every user in the safe operation of the nursing bed in accordance with these safety instructions, which must be supplied along with the unit.
- According to MPBetreibV 9 the user shall be notified of the location of this instruction manual.
- Call attention to the dangers of improper use of the unit, especially in regards to the electrical drives and side rails.
- The unit should only be operated by properly trained personnel.
- Make sure that representatives are also familiar with the operation of the nursing bed.
- During long term use it is recommended that the nursing bed be inspected for function and visible damage at regular intervals (Recommendation: yearly) (see chapter 6).
- Connect the bed directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.
- When attaching other equipment (e.g. compressors for positioning systems) make sure that they are mounted and can function in a secure and safe manner.
- Do not put multiple outlet electrical extensions under the nursing bed. Leaking liquids can cause a fire hazard.

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■ Pay special attention to:

Securely locating of all wiring, cables, tubes, etc.

- ▶ See page 10 : Intended use
- ▶ See page 22 : Safety note for the electric drive system

For further information please contact the manufacturer of the accessories or BURMEIER directly.

■ Make sure that your personnel follow the safety instructions.

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#### SAFETY INSTRUCTIONS FOR USERS

- Let the operator instruct you in the safe use of the nursing bed.
- Before each use make sure that the nursing bed is in proper and faultless condition.
- Take care that no obstacles, such as furniture or slanted ceilings interfere with the adjustment functions.
- Pay attention that when using additional electrical components, such as patient lifts, reading lights or compressors for positioning systems that their electrical cords do not get entangled or damaged by the moving parts of the bed.
- Connect the bed directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.
- Make sure that when attaching other appliances (e.g. compressors for positioning systems) that a safe mounting and function is guaranteed.
- Do not put multiple electrical outlets under the nursing bed. Leaking liquids can be a fire hazard.
- Pay special attention to securely locating all wiring, cables, tubes, etc.

For further information please contact the manufacturer of the accessories or BURMEIER directly.

- Take the nursing bed out of operation if damage or a malfunction is suspected:
- ▶ Immediately unplug the electrical plug from the wall outlet;
- ► Clearly mark the nursing bed as "DEFECTIVE";
- ▶ Promptly notify the appropriate operator.

A checklist for assessing the proper condition is located in Chapter 6.

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- Position the power cable in such a way that during normal operation of the nursing bed the cable will not be stretched, driven over or interfere with moving parts of the bed.
- Before any relocation of the nursing bed it is imperative that the power cable is pulled from the wall socket and be secured in the power cable holder.
- Check the power cable regularly for damage (abrasions, exposed wires, kinks, pressure marks, etc.) in particular:
- ▶ after every larger mechanical strain (e.g.: rolling over the power cable with the nursing bed or with an equipment cart; after a strong pulling or bending load caused by the bed rolling away with the power cable still attached to the wall outlet);
- ▶ after relocating and before plugging in the power cable.
- ▶ during prolonged operating by the user.
- Regularly inspect the cord grip of the power cable. It must be tightly screwed in.
- Do not use the area under the nursing bed as a "parking spot" for any utensils.
- Lower the reclining surface to its lowest position when leaving the patient unattended. This reduces the risk of injury to the patient getting in and out of the nursing bed.
- Keep the hand switch safe from accidentally falling during non-use (Hang it on the hook). Take care that the the cable is not damaged by the moving parts of the nursing bed.
- To protect the patient and especially children from unintentionally adjusting the hand switch, place it beyond their reach (e.g. at the foot end) or lock the hand switch when:
- ▶ the patient is not able to control the bed safely or is unable to get out of dangerous positions without help;
- ▶ the patient could be endangered, because of unintentional adjustments to the drive motor;
- the side rails are up: (There is a possible danger that limbs could be squeezed when adjusting the back and thigh rest);
- ▶ unsupervised children are in the room with the nursing bed.
- Adjustment functions can only be performed by a trained person!

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## PRODUCT DESCRIPTION NORMAL OPERATING PROCEDURE

- The nursing bed DALI II, hereafter called bed, serves as a comfortable solution for positioning and easing the care of infirm, frail people in senior or nursing homes. Additionally it was developed as a comfortable solution for in-home care of infirm, handicapped or frail people. It will support this care.
- The DALI II care bed is designed for use in hospitals only in rooms of application group 0.
- This bed can be used for care under the instructions of a physician and can serve in diagnosing, treatment and observation of a patient. It is therefore equipped with a locking mechanism for the electrical adjustment functions.
- This bed has no special connector possibilities for a potential equalization. Consider this before combining with additional electrically operated (medical) devices.

Further references about additionally applicable precautions can be found:

- ▶ In the instructions for the additional electrically operated devices (e.g. compressors for positioning systems, infusion pumps, nutrition probe ...)
- ▶ In the Standard DIN EN 60601-1-1:2002 (Safety of medical electrical systems)
- ▶ In the Standard VDE 0107:1994 (High current systems in hospitals)
- The maximum long-term load for this bed is 185 kg (patient and accessories). Under this load, the bed can be operated without restriction.
- The permitted weight of the patient depends on the total weight of the accessories attached at the same time (e.g. respirators, infusions,...)

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#### Example:

Weight of accessories (with mattress):	Max. permitted weight of patient				
10 kg	175 kg				
40 kg	145 kg				

- ▶ This bed is not appropriate for patients less than 150 cm tall. Observe, especially with patients in poor health condition, the safety tips in chapter 4.
- This bed may only be operated by trained personnel.

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- This bed is useful for multiple re-use. Bear in mind the necessary requirements:
- ▶ See chapter 5;
- ▶ See chapter 6.
- While the patient also is lying in bed, it is allowed to drive the bed within the building.

This bed may only be operated under the operating instructions laid out in this instruction manual.

Any other use is deemed to be inappropriate.

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#### STRUCTURAL DESIGN

The bed is shipped disassembled on a storage aid. It can be transported without problems even in apartment homes. It consists of two end boards (head- and foot section); a reclining surface frame which is divided in the middle; four side rail beams and a trapeze pole with trapeze handle. The bed stands on four steerable, individually-locking casters.

#### Reclining surface frame

The reclining surface frame has four sections: a movable backrest, a firm middle section, a moveable thigh rest and lower leg rest. The back- and thigh rest can be adjusted with the aid of electric motors. The horizontal height of the reclining surface can be adjusted and inclined to a low foot end position. All adjustment functions are controlled by a hand switch.

#### Side Rails

The bed has side rail beams on both sides which can be raised as a barrier or lowered when not needed. This safeguards the patient from accidentally falling out of the bed.

#### Electric adjustment system

The electrical adjustment system of this bed is error protected, flame-retardant (V0) and consists of:

■ the central control unit. For the safety of patient and user, safety low-voltage of 24 volt is produced in the control unit by means of a transformer. All drive motors and the hand switch are hooked up to the central control unit by plug-in connections and work with the 24 volt safety low-voltage;

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#### ■ a hand switch with a sturdy hook.

- The user can lock the adjustment functions of the hand switch if the poor clinical condition of the patient warrants is.
- two drive motors with integrated control for the horizontal height .
- $\blacksquare$  one drive motor with integrated control for the thigh rest .
- one drive motor with integrated control for the backrest.

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#### **DESCRIPTION OF MATERIALS**

The bed is mostly constructed out of steel sections, their surfaces have been polyester-powder-coated, or with a metallic coating of zinc and chrome. The end boards, the side rail beams and the spring wood slats consist of wood and derived timber products whose surface has been sealed.

All surfaces are harmless to skin contact.

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# ASSEMBLY AND START UP PROCEDURE

The bed ships dissasembled mounted on a storage aid.

The assembly is taken care of at the destination by qualified personnel supplied by the operator.

Remove all packaging material and cable ties.

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
#### **ASSEMBLY**

#### THE END BOARDS

The two end boards [3]+[10] must not be interchanged! The adjustment to a lower foot position can otherwise not be achieved.


Look for the decals marking the head and foot sections. They are located in the middle on the crosspipe, in the area of the drive motor mount and in the middle of the the crosspipe on the reclining surface frame [12].

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

#### ASSEMBLY OF THE RECLINING SURFACE


- Framove the side rail beams [7] and the trapeze pole [2] from the storage aid and put them aside for now.
- Remove the two halves of the reclining surface frame [12] from the storage aid.
- Tut the head side half of the reclining surface frame vertically on the floor. The two trapeze pole sockets [16] point downward while the 2 pairs of knurled screws are pointing upward.
- Now take the foot side half of the reclining surface frame and lift it over the head side half of the reclining surface frame. Now insert the halves of the reclining surface frame into each other.

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- Tighten the 4 knurled screws. Do not use pliers to tighten the screws.
- © Connect the lift rod of the drive motors with the entry of the head side of the reclining surface frame. For that, insert the cotter pins through the receiver and the lift rod and secure it by turning over the clamping yoke.
- Place the combined reclining surface frame flat on the floor.

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

#### ASSEMBLY OF THE END BOARDS

- Join the end board, head section [10] with the reclining surface frame [12] together. Ensure that the decals correspond!

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Lift the reclining surface frame at the head end and push the two support brackets from the end board, head section [10] into the tubes of the reclining surface frame [12] to the limit.

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There should only be a 5 mm gap between the reclining surface frame and the corner post of the end board.

- Tighten the 4 knurled screws tight. Do not use pliers to tighten the screws.
- Repeat the process with the end board, foot section [3].

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
#### **ASSEMBLY OF THE SIDE RAIL BEAMS**

At the end boards [3]+[10] you will find on the left side and right side a guide bar [17] for the side rail beams [7]. In them runs a plastic slide with 4 metal pins. The metal pins hook into the oblong holes on the face of the 4 side rail beams. The side rail beams can be raised or lowered, as needed.

- To mount the side rail beams the reclining surface frame has to be raised to mid-level (Electrical Connections see page 17).
- When first delivered there are two plastic slides in a guide bar.

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Start with the end board, head section:

- Toosen the knurled screws [a] on the bottom end of the guide bar [17] until you can remove one of the plastic slides toward the bottom.

  Leave the other plastic slide in the guide bar.
- In order to prevent the plastic slide from falling out of the guide bar, reinsert and tighten the knurled screws.

The plastic slide you have just removed must now be inserted into the second guide bar of the end board, head section.

- Loosen the knurled screws [a] at the bottom end of the guide bar [17] and insert the rounded end of the plastic slide into the guide bar.
- Reinsert and tighten the knurled screws.
- Pull both plastic slides on the head section upward until they click into place.

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#### Metal side rail

Now place the thin side rail bar [7] onto the top pair of metal plugs on the head board, and let the other end of the side rail bar rest on the floor.

The side rail beams are rounded on the top. The rounded side must always be uppermost.

Now place the thick side rail bar onto the bottom pair of metal plugs and let the other end of the side rail bar rest on the floor.

Please refer to the photograph.

Repeat the same procedure on the other side of the end board, head section.

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
#### Wooden side rail

Now place the side rail beam [7] onto the top pair of metal plugs on the head board, and let the other end of the side rail beam rest on the floor.

The side rail beams are rounded on the top. The rounded side must always be uppermost.

- Now place the second side rail beam onto the top pair of metal plugs and let the other end of the side rail beam rest on the floor.
- Repeat the same procedure on the other side of the end board, head section.

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
Now move to the end board, foot section:

Remove both plastic slides toward the bottom out of the guide bar.

### The rounded end of the plastic slide must face upward.

- Guide the plastic slide in the guide bar of the end board, foot section, upwards until it clicks into place.
- In order to prevent the plastic slides from falling out, reinsert and tighten the knurled screw [a].
- Repeat this process on the other side of the foot section.

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Check the distance between the beam and the guide bar. It should not be greater than 10 mm.

In case the distance is larger the reclining surface frame must again be pushed into place against the end boards.

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#### **ELECTRICAL CONNECTION**

The central control unit supplies the drive motors, through a spiral cable, with the needed voltage. The spiral cables from the drive motors must be inserted into the correct socket in the central control unit. Before connecting the cables all packaging materials must be removed.

Place the drive motor cables and the hand switch cable under the reclining surface frame.

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The spiral cables leading to each drive motor are marked in the following way at the central control unit:


- 1 Backrest motor
- 2 Thigh rest motor
- 3 End board / Head motor
- 4 End board / Foot motor
- **HB** Hand switch

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Take special care in assigning the drive motors for the height adjustments. These must not be interchanged. Spiral cable number 3 leads to the drive motor for the end board (head section). Spiral cable number 4 leads to the drive motor for the end board (foot section). Only one lower foot position should be possible.

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Make sure that no cable is damaged, that there is no loop in the cables and that the cables are not crushed by moving parts.

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When exchanging single electrical components, always make sure that their plugs have been plugged back into the central control unit until they lock and the plug covering strips have been re-affixed. This will ensure a tight fit and a proper operation.

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#### STRAIN RELIEF OF THE POWER CABLE

Underneath the reclining surface frame is the strain relief for the power cable.

Make sure that the power cable is hooked tightly into the strain relief.

Special care is to be taken when positioning the power cable. Make sure that the cable is not damaged and there is no loop in it! When moving the bed around, take care that the castors do not run over the power cable! When adjusting the bed the power cable should not get caught between moving parts!

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#### **LOCATION REQUIREMENTS**

- There must be enough space available to be able to adjust the bed. No furniture, window sills and the like should be in the way.
- The space under the bed must remain clear.
- Check before use on parquet flooring, that the castors do not cause discoloration of the varnish. Use on ceramic tiles, carpet, laminate and linoleum is harmless.

BURMEIER is not liable for any damage to floors due to daily operation.

- A properly installed 230 volt power outlet must be available, near the bed.
- Take care when attaching further accessory equipment (e.g. compressors for positioning systems, etc.) that their attachment is safe and their functionality is not impaired. Especially important is hereby also the safe placement of all movable connecting cables, tubes, etc. If you have any further questions or need clarification contact the manufacturer of the accessory equipment or BURMEIER.

Observe the following points when setting up the bed to minimize possible fire hazard. Instruct the user and the authorized personnel in these points!

- If possible use only fire-retardant mattresses and bed sheets.
- Only use appropriate firm mattresses according to DIN 13014. This mattress should be inflammable according to DIN 597 Section 1 and 2.
- Avoid smoking in bed because of the danger of fire.
- Use only technically faultless appliances (e.g. heating blankets) and other electronic devices (e.g. lamps, radios).

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- Take care that these devices are only used for the purpose intended and that they cannot accidentally get on top of or underneath the bed sheets (Fire hazard because of overheating).
- Absolutely avoid the use of extension cords or multi outlet extension cords under the bed (Leaking liquids can cause a fire hazard).
- Connect the bed directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.

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#### REACHING OPERATIONAL READINESS

- After successful assembly, complete an inspection according to chapter 6.
- Clean and disinfect the bed before first time use and before every further use according to chapter 5.

After successfully completing and heeding all steps on page 13 through 19 the bed is now operationally ready.

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
#### DISASSEMBLY

- Remove the trapeze pole.
- The Lower the bed into the middle horizontal position.
- Lock all four castors.
- Remove all side rail beams consecutively.
- Insert the plastic slide back into the guide bar. Screw the knurled screws tightly back into the guide bars.
- The Lower the bed into the lowest horizontal position.
- Pull the power plug out of the socket.
- Unplug the transformer cable from the connecting socket.
- Remove the plugs and strain relief from the drive motors at the head and foot ends.

### To disassemble the reclining surface frame proceed as follows:

- Lean the reclining surface frame upright on a wall. Head section downward.
- Remove the cotter pin from the lifting rod of the drive motor for the backrest.
- Loosen the knurled screws on both joints of the reclining surface frame.
- Pull the frame apart.
- Screw the loosened knurled screws back into the reclining surface frame so they do not get lost.

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#### STORAGE AID

The storage aid joins both headboards and takes in both halves of the reclining surface frame. Additionally there are holders for the side rail beams and the trapeze.

Screw both pieces of the storage aid onto the receiver pin of one of the end boards.

At the same time all knurled screws in the storage aid must face in one direction.

The receptacle for the reclining surface frame must point upward; the basket for the side rail beams go inward and the receptacles for the trapeze point toward the outside.

Use the knurled screws from the reclining surface.

- Screw on the second headboard.
- Put the first half of the reclining surface frame (foot section) on the longer receptacle so that the foot section points downward. The mattress brackets [14] point outward.
- Now put the second half of the reclining surface frame (back rest) onto the receptacle, so that the head section points downward. The mattress brackets point inward.

Safely secure the lifting rods of the drive motors with cable ties or something similar.

- Tighten the knurled screws.
- Put the side rail beams into the basket.
- Place the trapeze pole in the appropriate
- The bed is now ready for transport or storage.

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## SPECIAL SAFETY INSTRUCTIONS FOR THE DRIVE UNIT

- Take care when adjusting the bed that no limbs of patients, users or any other persons, especially children, can be jammed or injured by getting under the movable parts or the bed frame.
- For the safety of the patient, lock the hand switch operation of adjustments to the back and thigh rests to avoid unwanted adjustments when the side rail beams are up. (Danger of getting limbs trapped when adjusting the back- and thigh rests).
- When using accessories with the bed it is necessary to: make sure that the arrangement of these accessories will not cause any places that might squeeze or hurt the patient. If this is not the case the operator must safely disable the adjustment function of the back- and thigh rests.
- ▶ To do so place the hand switch out of the patient's reach (e.g. foot section) or lock the adjustment functions on the hand switch.
- Make sure that the power and hand switch cable cannot be jammed or damaged in other ways.
- Make sure that before relocating the bed, the power cable cannot be stretched, rolled over or damaged in any other way.

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- If the maximum continuous run time is not heeded (e.g. by constantly "playing" with the hand switch) a thermal safety component switches off the power for safety reasons.
- The adjustable range of every function is electrically/mechanically possible only within a predetermined range.
- Electrical disturbances affecting other nearby electronic appliances (e.g. crackling on the radio) cannot be excluded, despite adherence to all prescribed regulations. In such rare circumstances increase the distance between the appliances, do not use the same power socket or switch the affected appliance off for a while.

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HAND SWITCH

The electrical adjustments can be made using the hand switch [8].

The adjustable range of every function is electrically/mechanically possible only within a predetermined range.

Locking mechanisms have been built into the hand switch for safety reasons. For the safety of the patient, adjustment functions can be locked at the hand switch, if, from a clinical point of view and with the approval of the attending physician, the condition of the patient warrants

- The hand switch, with its elastic hook can be hooked onto the bed.
- The spiral cable allows for greater mobility.
- The hand switch is water-tight and washable (IP X4).

#### **OPERATING THE HAND SWITCH**

- Train the patient in the functions of the hand switch.
- The drive motors will run as long as the button is held down.
- The adjustment functions, except for the lowered foot position, can be used in both directions.
- Basically the functions of the buttons are as shown:
  - raise
  - lower

#### Adjusting the backrest

Push this button if you want to adjust the backrest angle.

#### Adjusting the reclining surface height

Push this button if you want to adjust the height of the reclining surface.

If the bed is in the lowered foot position, it will automatically be adjusted to horizontal level, when it is readjusted to the lowest or highest position.

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#### Adjusting the thigh rest

Push this button to adjust the thigh rest angle.

#### Adjusting to the lowered foot position

Push this button if you want to adjust to the lowered foot position. With this button ▼ only the function "Down" is possible.

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■ Unlock the foot brake of both castors on the head or foot section to avoid damaging the floor.

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LOCKING FUNCTION

### Only the operator is authorized to use the locking function!

In case the clinical condition of the patient is so critical that any adjustment of the beds by means of hand switch causes danger, the operator has to disable the bed immediately. The nursing bed remains in the position it had at the time of disconnecting.

Do not turn the locking key with force over the stop of the lock! The lock or the whole hand switch can be damaged.

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The locking key is delivered attached to the bed with a cable tie on the hand switch. The locking key is not intended for use by the patient. It must be removed from the hand switch. The operator or a person authorized by the physician should keep the **locking key**.

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■ Turn the respective lock on the hand switch with the locking key clockwise in the locking position.

The color of the respective display will change from green to yellow.

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Lock is vertical

Color of the display is: green

Buttons can be pushed ("Clicking sound")

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
#### **Drive disabled**

Lock has been turned clockwise by about 15°

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#### **CASTORS**

The bed stands on four steerable, individually-locking castors [11].

- Braking: Push the foot lever down.
- **Driving**: Push the foot lever toward the corner pole.
- Relocate the bed only when it is positioned in the lowest position.
- Before relocating the bed, ensure that the transformer unit is securely positioned on the bed so that it cannot fall off.

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#### **POWER CABLE HOLDER**

The power cable holder is located sideways on the reclining surface frame near the foot section.

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
■ The maximum load of the trapeze pole is 75 kg at the front end.

The trapeze pole attached to the bed makes it easier for the patient to get into and out of bed.

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#### Anbringen

Den Aufrichter [2] in die Hülse stecken. Der Metallstift [b] muss sich in der Aussparung [c] der Hülse befinden. In both corners of the head section of the reclining surface frame are two trapeze pole sockets [16]. The trapeze pole socket has a notch at the top [c], which together with the metal pin [b] restricts the swivel range of the trapeze pole. The trapeze pole should be installed on the side where the patient gets in and out of bed.


#### **Attaching**

Insert the trapeze pole [2] into the socket. The metal pin [b] must be in the notch [c] of the socket.

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### Removal

Pull the trapeze pole straight up out of the socket.

#### Trapeze pole swivel range

Do not swivel the trapeze pole outside of the bed. There is a danger that through the pull, the bed will tip over.

The metal pin must always be in the notch of the socket.

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#### **TRAPEZE**

A trapeze (grab handle) can be attached to the trapeze pole. The patient can hold on to the trapeze to make it easier to change to a different position.

Regularly check the trapeze at the strap for any damage (see Chapter 6). A damaged trapeze or strap must immediately be exchanged.

#### **Durability**

There is a date clock on the trapeze. With normal use the trapeze will last for at least five years. Afterwards a sight check is to be performed periodically to determine if further use is permissible.

#### Adjustable range of the trapeze

The height of the trapeze can be adjusted by means of the strap to accomodate a range of about 55 cm to 70 cm (measuring from the top of the mattress).

- Push the loop of the trapeze [1] over the first stud of the patient grab handle pole.
- Test by tightly pulling the trapeze downward that it is secure.

**Useful Tip:** The maximum lifting load of the trapeze pole is 75 kg at the front end.

- The height of the trapeze can be adjusted by means of the strap.
- Make sure the the strap is fed correctly through the buckle.
- Make sure that the end of the strap is at least 3 cm above the buckle.

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#### **LOWER LEG REST (OPTIONAL)**

When the thigh rest is fixed, the lower leg rest [5] can be individually adjusted. For this reason there are two adjusting ratchets, which will lock the lower leg rest into position.

#### Raising

Raise the lower leg rest on the foot end - not with the mattress brackets - up to the desired level. The lower leg rest will automatically snap into place.

#### Lowering

- First lift the lower leg rest all the way to the stopping point.
- Slowly lower the lower leg rest.

There is a potential for injury if the lower leg rest falls straight down.

If you lower the thigh rest, the lower leg rest will be lowered automatically.

#### SIDE RAILS

Side rails conveniently prevent the patient from falling out of bed accidentally. They do not, however, prevent a patient from leaving the bed.

#### Engaging

- Pull the side rail beams [3] on one end after the other until until they lock in its highest position. You should not be able to move them up or down.
- Theck that the side rail beams are locked into place by pushing on them from above.

#### Disengaging

- Lift the side rail beams slightly.
- Push the upper release lever [4] and slowly lower the side rail beams.

### Special safety notes for use of side rail beams:

- Only use safely lockable, technically perfect, undamaged side rail beams with proper gap dimensions
- Make certain that only suitable Original BURMEIER side rail beams are used, which have been approved for use with the bed model you are using or that are already integrated in the bed.

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- Before using side rail beams evaluate and take into account the clinical condition and any special problems of the patient
- ▶ If, for example, the patient is very bewildered or very anxious, forgo the side rails and use alternative safety measures such as a restraining belt or restraint sheets.
- ▶ For especially small, frail patients, additional protective measures for reducing the distance between bars on the side rails may be necessary. In these cases use, for example side rail foam covers (accessories, see chapter 8), retaining belts, etc.

This will effectively ensure patient safety by reducing the risk of becoming trapped or slipping through.

- Only use appropriate, firm mattresses accord to DIN 13014 with a height of 10 to a maximum of 18 cm. Further, this mattress should be fire resistant according to DIN 597 Section 1 and 2.
- If elevated patient positioning systems are employed, such as antidecubitus mattresses, an effective side rail height of at least 22 cm above the unweighted mattresses must be ensured. If this requirement is not met, the operator must take suitable measures to safeguard the patient. This can be done by using an attachable bar for side rails (Accessory).
- The mechanical parts on the bed frame and the side rails which serve to protect the patient must be regularly examined for damage. This should not only happen before the installation of the side rail but also during use in shorter intervals (at least before every new assignment).
- To safeguard the patient against unintentional adjustments, always lock the electrical adjustment of the back and thigh rests on electrical beds when the side rails are raised. Otherwise there is a risk of trapping limbs when the back and thigh rest are adjusted.
- Before lowering the side rail beams make sure that the patient's limbs are clear of the side rail beams. This way you can avoid injuring the patient.

If these safety warnings are not observed there is a danger for the patient of being injured by the side rails, of falling out of bed because of too large a distance between the bars or because of damage, the use of unsuitable side rails or the incomplete locking of the side rail.

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
## EMERGENCY LOWERING OF THE BACK REST

The upright back rest can be lowered manually in case of a power outage or failure of the power supply or the electrical drive systems.


In this case two users are absolutely necessary!

Disregarding these safety guidelines and instructions may cause the backrest to fall out of control, which could lead to serious injuries for the user and the patient!

Emergency lowering may only be carried out in the case of extreme emergency and by users who have a complete command of the procedures below.

We strongly advise you to practice emergency lowering of the backrest under normal conditions so that in case of an emergency you will be able to react quickly and properly.

- Before the emergency lowering, ease the load on the backrest.
- The first user raised the back rest slightly by the outside edge of the head section and holds it in this position.

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To convert the bed back to its normal operating state: Swing up the connecting mount and secure the spring pin, then turn over the bent bracket. Kapitel 5

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# CLEANING AND DISINFECTING

#### **GENERAL**

Cleaning is the most important requirement for a successful chemical disinfection.

Generally a routine cleaning of the bed, if used by the same patient, is hygienically sufficient.

Disinfection of the bed frame is only necessary in case of visible contamination with infectious or potentially infectious material (blood, stool, pus) or by presence of an infectious disease under the direction of a physician.

The bed must be cleaned and washed down with disinfectant in case of a patient change!

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### Before beginning the cleaning, please note:

- Unplug the power cable from the wall socket and store it in such a way that it does not come into excessive contact with water or other cleaning solutions (store in plastic bag).
- Ensure that all plugs on the transformer unit and the drive motors are plugged in properly.
- None of the electrical components must show any visible damage; otherwise water or cleaning solution could penetrate. This could lead to malfunction or damage to the electrical components.
- Before operating the bed again, make sure that there is no residual moisture on the electrical contacts, by drying or blowing on the power plug.
- The electrical components must not be cleaned with a water jet, a high pressure cleaner or other similar device!
  Clean only with a moist towel.
- If you suspect that water or any other form of moisture has penetrated the electrical components, pull the power plug out of the wall socket immediately. If already disconnected from the power supply, make sure it is not plugged in again. Clearly mark the nursing bed as "Defective" and take it immediately out of service. Report this to the operator promptly.
- Failure to follow these directions could lead to major damage of the bed and its electrical components and lead to subsequent malfunctions.

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### CLEANING AND DISINFECTION

- Remove the bed sheets and give them to the laundry.
- Clean all surfaces including the slatted boards and the reclining surface made of plastic inserts or wire mesh base with a mild and environmentally safe cleaning solution. This also applies to the hand switch.
- If the bed has visible contamination e.g. infectious or potentially infectious material has polluted it, it must be disinfected by wiping after cleaning. Use approved disinfectants for each surface according to the products listed by the DGHM (German Organization for Hygiene and Microbiology). This also applies to all beds from patients with illnesses subject to report according to §6 of the infection control law (German abbreviation: IfSG, Protection against Infection Act), colonization or infection with treatment resistant pathogens (e.g. MRSA, VRE) and all beds from the intensive care and infection control units. Use only the DGHM approved disinfectant in the approved concentration.
- Disinfecting the castors is only necessary when visibly contaminated by infectious or potentially infectious material.

**Notice:** Continuous disinfection is only necessary for patients with a treatment resistant pathogen (e.g. MRSA) within the hospital.

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## INSTRUCTING THE USER AND QUALIFIED PERSONNEL

To guarantee the correct procedure for cleaning and disinfection is followed, we recommend to correctly instruct the handler and qualified personnel. Thereby communicating that the following points need to be observed:

- The clean bed is to be transported to the patient's home in such a way that it does not become polluted or contaminated.
- When disassembling the bed we recommend cleaning and disinfecting it right away. The qualified personnel should be familiar with the steps necessary to thoroughly clean and disinfect (the operator will provide the individual steps necessary). It is necessary to use only the DGHM approved disinfectant in the approved concentration. The disinfectant must be appropriate for the surface being cleaned.
- The qualified personnel must be equipped with water-tight one-time-use aprons and gloves.
- Only fresh and clean towels must be used to clean and disinfect, these must then be given to the laundry.

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■ After cleaning the bed, the qualified personnel must also disinfect their hands before going to the next patient. An appropriate disinfectant dispenser should be part of the qualified personnel's gear.

The immediate cleaning of the bed while disassembling has the benefit of not transporting "dirty" beds or bed parts along with clean ones in the same truck. A transfer of potentially infectious spores which may be clinging to the used bed may be avoided. A transfer of spores of a hospital-acquired infection will be avoided with strict adherence to these recommendations.

If the beds are to be stored for a while, they should be covered to keep them free of dust, accidental pollution and contamination.

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## CLEANING SUPPLIES AND DISINFECTANTS

Make use of these recommendations to keep the bed serviceable for a long time.

- The surfaces must be intact; any damage should be repaired immediately.
- We recommend a (damp) wash down. When choosing the cleaning solution make sure it is mild (gentle to skin and surfaces) und environmentally friendly. In general a regular household cleaner can be used.
- For cleaning the coated metal parts use a wet cloth with regular household cleaners.
- The disinfectant used for the wash down should be one approved for the purpose and listed by the DGHM (in the recommended dosage). Disinfectants and concentration listed by the Robert-Koch-Institute (the former Bundesgesundheitsamtes) should only be used in case of an epidemic according to the infection safety law and only on the orders of the appropriate health department.
- Despite the excellent mechanical resistance of the bed, scratches, markings, etc. which penetrate the coating should be resealed using a suitable medium to prevent the penetration of moisture. Contact BURMEIER or your specialized dealer with any questions.

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■ Do not use scouring agents, stainless steel cleaners, abrasive cleaning agents or scouring pads. These substances can damage the surface.

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Notice:

Surface disinfectants based on aldehyde generally have the advantage of working on a large array of pathogens, have a relatively low protein error and are environmentally friendly. Its main disadvantage is the allergy-causing and irritation potential. Compounds on the basis of glucoprotamine do not have this disadvantage and are very effective, though mostly more expensive. Disinfectants on the basis of chlorine splitting combinations, with prolonged contact or in too high doses, can be corrosive to metals, plastics, rubber and other materials. In addition to this, the compounds have a high protein error, are irritants to the mucous membranes and are not environmentally friendly.

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#### **CONTACT WITH DISINFECTANTS**

- Follow the guidelines for each product and manufacturer. Observe the exact dosage! We recommend automatic dosage apparatus.
- Always mix the solution with cold water, to avoid vapors that aggravate the mucous membranes.
- Wear gloves to avoid direct skin contact.
- Do not store the mixed disinfectant in open containers with the cleaning rag inside. Close the container!
- Use closeable bottles with dosage pumps for wetting the cleaning rag.
- After finishing the surface disinfection ventilate the room.
- Use a wipe down disinfecting method, do not use a spray. Using the disinfectant as a spray causes the disinfectant to turn into a mist which could be inhaled.
- The wipe-down disinfection plays a major
- Do not use alcohol for larger areas.

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
operators of nursing beds are obligated to keep medical products in a safe state throughout their entire service life. This also includes regularly carrying out expert maintenance as well as safety checks.

If a bed is purchased for domestic use (not for commercial use), the regular safety inspection is not compulsory but is still recommended by the manufacturer.

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#### **Notice for Operators**

This bed has been constructed and built to last safely for a long time. With appropriate operation and use, this bed has a life expectancy of 2 to 8 years. The life expectancy is based on the conditions and frequency of use.

Through repeated transport, assembly and disassembly, improper operation as well as long time use, it is possible that damage, defects and deterioration can occur. These deficiencies can lead to danger if they are not recognized and quickly fixed.

Because of this there are legal bases for the regular enforcement of tests to ensure the safe condition of this medical product.

According to § 4 of the Medizinprodukte-Betreiberverordnung operators are obliged to maintain the product. Because of this the following regular inspection and function controls are to be performed by the operator as well as the user.

Point out these necessary tests to the user.

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#### BY THE USER

Along with the regular substantial tests by technically qualified personnel, the normal user (nurses, supporting relatives, etc.) must also perform a sight and function check at shorter, regular intervals, as well as before a new patient.

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#### Recommendation:

All electrical and mechanical components should be tested monthly. The main power cable as well as the hand switch cable should be checked after every mechanical load and after every relocation.

For this purpose use the checklist on page 36.

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### Checklist: Test protocol by the user

Test		ОК	Not OK	Defect description
Visual check of electrical components				
Hand switch and cable	Damage, Cable Routing			
Power cord	Damage, Cable Routing			
Hand switch	Damage, Foil			
Transformer unit	Damage, Cable Routing			
Visual check of mechanical componen	ts		•	
Trapeze pole/Socket	Damage, Deformation			
Bed frame	Damage, Deformation			
Spring wood slats	Damage, Splinters			
Bed box	Damage, Splinters			
Reclining surface frame	Damage, Deformation			
Side rail	Damage, Splinters			
Function test of electrical components				
Hand switch	Functiontest, Blocking function			
Function test of mechanical componer	nts			
Emergency lowering of backrest	Test according to instruction manual			
Castors	Free movement an locking			
Knurled screws	Tight fit			
Side rail	Safe locking/unlocking			
Lower leg rest	Latching			
Accessories (e.g. trapeze pole, trapeze)	Secure, Damage			
Signature of the tester:	Result of the test:			Date:

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#### BY THE OPERATOR

The operator of the nursing bed according to MPBetreibV § 4 is obligated with every new assembly, every maintenance and throughout the operation, to perform regular inspections, to ensure the safe condition of the nursing bed.

These tests are to be repeated during regular maintenance and depending on the use according to MPBetreibV § 4 and from the Directive of the German Employers Liability Insurance Association according to BGV A3.

- Observe this sequence for the inspection according to DIN EN 62353:
- I. Visual inspection
- II. Electric Measurements
- III. Function Test
- ◆ The performance of the function test as well as the rating and documentation of the test results can, according to MPBetreibV § 4 only be performed by experts with the appropriate knowledge, training and experience as well as the tools necessary to perform it according to the rules.
- ◆ The performance of electrical measurements according to DIN EN 62353 can, with availability of proper measuring equipment, also be performed by an electro-technically trained person [according to BGV A3] with medical and appliance specific training.
- The assessment and documentation of the test results can only be performed by an expert with the corresponding knowledge, training and experience.

#### Test method:

- Output power test: direct or differential according to DIN EN 62353
- Output power test according to the test unit manufacturers instructions.

#### Limit:

■ Output power less than 0.1mA.

#### Test cycle:

We recommend a yearly visual inspection and performance check. If this checks was been ok, an electrical testing suffices every three years.

If the 2% error quota (see BGV A3: §5, table 1 B) is shown to have been observed, the test cycle for inspections and functional testing can be extended to a maximum of two years.

If damage or a function error is suspected take the slatted frame out of service immediately and disconnect it from the power until an exchange or repair can be made!

The test protocol draft on **pages 40 and 41** can be used.

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Customer / med. Facility / practice:							
Address:							
Carried out: ☐ Repeat inspection ☐ Inspection prior initial operation					ation		
□ □ Inspection following repairs / servicing							
Equipment type: ☐ Hospital bed ☐ Nursing bed Pr			Protection classification:				
Bed type: DALI II			Inventory Number:				
Location:			Serial Number:				
Grade level (IEC60601-2-52): □ 1 □ 2 □ 3 □ 4 □ 5							
Manufacturer: Burmeier GmbH & Co. KG			User-specific: none				
Testing equipment used (Type / Inventory Number):			1.				
MPG Classification: Class I			2.				
I. Visual Inspection				ок	Not OK	Defect Description	
Visual inspection of the electriacl components							
Stickers and type plates	Presence, legib	le					
Housing of transformer unit / Power adapter	Damage						
Motor housing and lifting rods of motors	Correct position	n, Damage					
Hand switch	Damage, Foil in	ntact					
Motor-, hand switch- and power cables	Damage, Routing						
Plug and plug covering strip of control unit	Available, Correct position						
Visual inspection of the mechanical components							
Stickers and type plates	Presence, legible						
Safe working load / patient weight	Presence, legible						
Trapeze pole / socket	Damage, Deformation						
Bed frame	Damage, Deformation						
Spring wood slats	Damage, Splinters						
Castors	Damage						
Reclining surface	Damage, Deformation						
Wooden boards	Damage, Splinters						
Welded seams	Split welded seams						
Side rail	Damage, Splinters						
Knurled screws	Tight fit						
Wearing parts such as joints	Damage						
II. Electric measurments according to DIN EN 62353:2008-08 (leakage current, direct)							
Proceed as follow:				ОК	Not OK	Defect Description	
Plug the power cable into the test socket on the measuring instrument.							
2. Connect the measuring instrument probe with an exposed conductive part of the bed frame (screws, etc.)							
3. Activate the drive motors by using the hand switch for the duration of the measurement.							
4. Start procedures with the measuring instrument							
	Limit Value	Value of first measurement	Current actual value				
Result: Bed SK II (Type B)	0,1 mA	mA	mA				


Test protocol of an inspection of electro-medical equipment according to DIN EN 62353: 2008-08 Page 1 of 2

Not OK **Defect Description** III. Performance check OK Performance check of the electrial components Motion limit switching of the motors Automatic cut-off Hand switch, blocking functions Test according to instructions for use Control unit and drive motors Abnormal Noise Development Control unit and drive motors Test according to instructions for use Strain relief for power cable Power cable firmly fastened Performance check of the mechanical components Joints and pivots Smooth operation Trapeze with belt strap Secure hold Castors Drive and brake Emergency lowering of the backrest Test according to instructions for use Side rail Safe locking / Secure hold / Safe unlocking Lower leg rest Latching Accessories (e. g. trapeze pole, trapeze) Secure, Damage, Suitability Result of the inspection: All values within permissible range: ☐ yes □ no Next inspection date: Inspection passed: □ yes □ no If inspection was not passed: ☐ Defective, do not use the bed! => Repair ☐ Defective, do not use the bed! => Take out of Service  $\hfill\square$  Bed does not meet the safety standards Test approval sticker applied: ☐ yes □ no Part of this inspection report: Comments: Signature: Inspected on: Inspector: Evaluated on: Operator / Expert: Signature:

Test protocol of an inspection of electro-medical equipment according to DIN EN 62353: 2008-08 Page 2 of 2

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#### **SPARE PARTS**

The corresponding replacement parts can be obtained from BURMEIER by specifying the item number, the order number and the serial number. The necessary details are found on the type plate under the cross tubing of the reclining surface frame.

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#### MANUFACTURER'S ADDRESS

In order to maintain operational safety and the right to claim under warranty, only original BURMEIER replacement parts may be used!

For ordering replacement parts, customer service requests and other questions, please contact:

**Burmeier GmbH & Co. KG** (A subsidiary of the Stiegelmeyer Group)

Pivitsheider Strasse 270 32791 Lage/Lippe, Germany Phone: ++49 52 32 / 98 41-0 Telefax: ++49 52 32 / 98 41-41

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# REPLACEMENT OF ELECTRIAL COMPONENTS

# Dangerous voltage, electric shock possible

- Before doing any work unplug the main power cable from the wall socket!
- Any work on the electrical equipment may only be performed by customer service, the motor manufacturer or qualified and authorized electricians under consideration of the applicable VDE rules and safety guidelines!

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#### **Exchange of the Hand switch**

- Unplug the power plug from the wall socket.
- Unplug the plug at the end of the spiral cable of the hand switch from the connector underneath the reclining surface.
- Lay the cable of the new hand switch so that it cannot be damaged by moving parts of the bed.
- Insert the plug at the end of the spiral cable of the new hand switch into the connector of the hand switch.

Take care that the O-ring on the plug is not damaged, otherwise the connection may not be moisture-proof.

Tarry out a function test of the electrical adjustments after the exchange.

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#### TROUBLESHOOTING

The following table on page 44 is a guide to rectifying common malfunctions:

If a function fault is detected during operation, which the user cannot rectify with the use of the trouble-shooting table aid, the specialized staff for maintenance and repair of the appropriate owner is to be called.

The user may **under no circumstances** try to eliminate defects at the electrical equipment!

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# Dangerous voltage, electric shock possible!

Any work on the electrical equipment may only be performed by service personnel, the motor manufacturer or qualified and authorized electricians as required by the applicable VDE rules and safety guidelines!

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#### Trouble-shooting table aid

Problem	Possible causes	Solution		
Hand switch / drive system not functioning	■ Main power cable is not plugged in	☞ Insert main power cable		
	■ No power supply to socket			
	■ Hand switch plug not properly connected	Test the plug connection		
	■ Hand switch or drive system defective	Inform your operator about necessary repairs		
	■ Functions are blocked on hand switch	Unlock functions (see page 24)		
Drives only operates for	■ Too much weight on the bed			
a short time when buttons are pressed	■ Bed is blocked by an obstacle			
Individual drives operate in one direction only	■ Defective hand switch, motor drives or control unit			
The wrong function is active when hand switch buttons are pressed	■ Internal motor plugs incorrectly connected (mixed up)			
Drives stop suddenly after prolonged adjust- ment times	■ Thermo switch in transformer (drive unit) is triggered on	Pay attention minimum pause (see page 22)		

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### **ACCESSORIES**

Only the Original BURMEIER accessories may be used. If other accessories are used, we can not be liable for any resulting injuries, defects and dangers.

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Make sure that the arrangement of accessories does not produce any crush or shearing zones for the patient when the back and leg rests are adjusted. If this cannot be ensured, the user must safely prevent the patient from adjusting the back- and leg rest.

▶ To do so, attach the hand switch out of reach (e.g. at the foot end of the bed), or lock hand switch adjustments.

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# MEASUREMENTS AND DIMENSIONS

#### Assembled bed with side rails:

Reclining surface External dimensions

Total weight, depending on model

Safe load

Patient weight, max.

#### Disassembled bed:

End boards

Reclining surface frame with motors

4 side rail bars, wooden 4 side rail bars, iron

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#### **SPECIAL FEATURES**

- Electrical height adjustment of the reclining surface from about 40 cm to 80 cm.
- Electrical adjustment of the backrest from 0° to about 70°.
- Electrical adjustment of the leg rest from 0° to about 35°.
- Electrical inclination of the reclining surface to the foot end from about 10°.
- Driven with four individually-locking castors.
- Reclining surface of 200 x 90 cm, divided into four areas; Outside dimensions are about 212 x 101 cm
- Side rails on both sides can be lowered.

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#### **ELECTRICAL DATA**

#### Control unit

(flange mounted to the backrest motor)

Type

Input Voltage

Max Amps

**Output Volts** 

**Output Amps** 

Operating time

Protection, splash proof

Classification, Protection class II, Type B, MPG classification Class I, not for use in explosive atmospheres.

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Storage temperature

Relative humidity

Noise level during motor adjustments

The ambient conditions stated below must be

Medical Products Directive

Repeat testing ... of medical electrical units

Technical aids for handicapped persons

Handicapped friendly construction

Risk analysis for medical products

Adjustable beds for handicapped persons

Safety for electro-medical equipment

Electromagnetic tolerances

Electrically operated hospital beds (reduced to applications for in home care)

Medical electrical equipment - Particular requirements for medical beds

Classified as class I active medical product (in accordance with MPG § 13).

Fulfills all additional safety requirements by the Highest State Authority from May 22, 2001 for

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### **DISPOSAL INSTRUCTIONS**

- This nursing bed is if electrically adjustable classified as a commercially used electric appliance (b2b) in accordance with WEEE-Guidelines 2002/96/EG (Electronics Law).
- Exchanged electronic components (drives, control unit, hand switch, etc.) of these nursing beds are considered electronic scrap and must be disposed of according to the WEEE-Guidelines.
- ◆ With nursing beds that came into use after August 13, 2005 (see designation "Serial" on the type label at the head of the bed), the operator is required to not dispose of the electrical components in the municipal waste area, but to send them directly to the manufacturer. BURMEIER and their service and sales partners will take these parts back.
- Our general terms and conditions apply for this redemption.
- No longer useable rechargeable batteries have to be disposed of according to regulations and do not belong in the house garbage.
- With all disposed components the owner must make sure that they are not infectious or contaminated.
- In the event that the nursing bed must be scrapped, the plastic and metal parts used must be disposed of separately, following the proper disposal procedure.
- Pneumatic springs are under high pressure! Before disposal they have to be depressurized. Please contact the manufacturer of the device on how to depressurize them (see type plate).
- Hydraulic units are made up of mostly metal and plastic. Before disposing according to manufacturer guidelines the oil has to be drained and disposed of properly.
- Contact your local community, disposal facility or our service department with any questions.

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We,

Burmeier GmbH & Co. KG Pivitsheider Strasse 270 32791 Lage / Lippe Germany

hereby declare that the product named below:

### Nursing bed DALI II

complies with the regulations of EC Guidelines 93/42/EWG for medical products and fulfills all additional safety requirements by the Highest State Authority from May 22, 2001 for care beds.

#### Applied Standards:

DIN EN 62353: 2008-08 Repeat testing of ... medical electrical units
 EN 12182:1999-11 Technical aids for handicapped people
 DIN 32977-1:1992-07 Handicapped fair constructions
 DIN EN 1970: 2005-10 Adjustable beds for handicapped people
 DIN EN 60601-1: 1996-03 Safety for electro-medical equipment
 EN 60601-1-2:2002-10 Electromagnetic tolerances
 IEC 60601-2-52: 2009 Medical electrical equipment

 (part 2-52: Particular requirements for the safety and essential performance of medical beds)

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Email

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No part of this manual may be reproduced without the prior written permission of the publisher!

All rights reserved.

Subject to technical changes!

 $\epsilon$ 

Stand: 2010.09.20 Updated: 2010.09.20

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Burmeier GmbH & Co. KG (Ein Unternehmen der Stiegelmeyer-Gruppe)